

2013 - 2014 Product Guide



Dansereau

RX SOLUCIONES DE IMAGENOLOGÍA S.A. de C.V.
Av. Moctezuma 3515-1B, Col. Ciudad del Sol. Zapopan, Jalisco. MÉXICO

VISITE Nuestro Sitio Web
y vea nuestra Línea de
Productos Completa



TEL (33) 3880-1262 CEL (33) 3499-7095 / 1300-5813

contacto@rxdigitalsoluciones.com
ventas@rxdigitalsoluciones.com

RXDigital

RX Soluciones de Imagenología

www.rxdigitalsoluciones.com

Dansereau Dental

Dansereau has over 50 years experience in manufacturing a full line of quality dental products and making them available directly to the dentist. Our Customer Support is exceptional and adds to the benefit of choosing Dansereau as your equipment supplier. Dansereau is committed to becoming your equipment choice for the life of your practice.

Warranty

By choosing Dansereau as your equipment supplier you will not only receive quality American made dental equipment, you are also entitled to our "Industry Leading 10 Year Warranty".

Support

Simple designs, quality parts and care to the details equal a superior product with an outstanding warranty.

Technology Support:

Dansereau has moved customer service into the 21st Century with a website that the dentist, staff or any support technician can use from a computer, tablet or Smart Phone, it lists all the parts necessary to support your Dansereau equipment as well as a full series of on-line videos for maintenance.

Advanced Video Support:

Dansereau has created a new video conferencing system that allows dentists, support staff and service technicians direct access to the factory via "Skype" or "Tango". You can use your Smart Phone or tablet to video conference with Dansereau factory personal for maintenance and service support.

Continuity

Dansereau has our Customer Service Department contact all our dentists every 90 days to verify the operation of the equipment is meeting the dentists standards and inquire if any maintenance parts are required.




NOTE:

In order that product improvements may be made at any time, specifications and other data are subject to change without notice or obligation to modify previously manufactured items. Misprints and errors in invoicing will not oblige Dansereau to honor a misprint or error in cost.

Pure American Reliability

Dansereau is a complete manufacturing facility with internal machining centers to fabricate their parts. They have designed all their equipment around the concept of reliability, ease of use and user friendly maintenance. This all being confirmed by our industry leading 10 Year Warranty.

- Internal Machining Centers Ensure Quality Parts



- Factory Assembly
- Guaranteed Reliability
- Ensure Quality



10 Year
Warranty

- Custom Fabricated Cushions with 100 Color Choices



Multiple Cushion Systems

- Memory Foam Cushions
- Memory Aseptic Cushions
- Aseptic Cushions - Seamless
- Removable / Reversible

- Dansereau **does not** outsource their cushion fabrication, all upholstery is done in factory, offering the dentist more choices.

INFINITY DENTAL CHAIR

- Ultra Thin Seat Back
- “**Memory Foam**” Cushions
- Dual Hands Controls Operation
- Auto Return & Adjustable Preset
- Drop Down Arms for EZ Access
- 800lb Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

10 Year
Warranty



AM3 CONTEMPORARY DENTAL CHAIR

- Ultra Thin Seat Back
- “**Memory Foam**” Cushions
- Dual Hands Controls Operation
- Auto Return & Adjustable Preset
- Drop Down Arms for EZ Access
- 800lb Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

10 Year
Warranty



AM3 TRADITIONAL DENTAL CHAIR

- Ultra Thin Seat Back
- Aseptic Seamless Cushions
- Dual Hands Controls Operation
- Auto Return & Adjustable Preset
- Drop Down Arms for EZ Access
- 800lb Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

10 Year
Warranty



ACCESSORIES

Chair Foot Control -

Chair Swivel -

Articulating Headrest -

AM SURGERY DENTAL CHAIR

- Ultra Thin Seat Back
- Increased Chair Height - 27"
- Dual Hands Controls Operation
- Auto Return & Adjustable Preset
- Drop Down Arm Boards for EZ Access
- 800lb Lift Capacity
- Over 150 Color Choices

10 Year
Warranty

AM PEDO DENTAL CHAIR

- Ultra Thin Seat Back
- Aseptic Seamless Cushions
- Dual Hands Controls Operation
- Auto Return & Adjustable Preset
- Drop Down Arms for EZ Access
- 800lb Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

10 Year
Warranty

CALIFORNIA DENTAL CHAIR

- Thin Seat Back
- Removable/Reversible Cushions
- Dual Hands Controls Operation
- Automatic Return Included
- 800lb Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

10 Year
Warranty

ACCESSORIES

Chair Foot Control -

Chair Swivel -

Articulating Headrest -



America Value Operator

10 Year
Warranty

AMERICA VALUE OPERATORY

- AM3 Dental Chair
- 3 HPC Unit and 3 Way Syringe
- Telescoping Assistants Package
- Diamond Dental Light
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY



COMPLETE LIST OF ACCESSORIES ON PAGE 17



Universal Value Operator

10 Year
Warranty

UNIVERSAL VALUE OPERATORY

- AM3 Dental Chair
- Radius 3 HPC Cart Unit and 3 Way Syringe
- Radius Telescoping Assistants Package
- U Cart with Adjustable Height
- Laminated Worksurface Area
- Diamond Dental Light
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY



COMPLETE LIST OF ACCESSORIES ON PAGE 17



Clinic Operator



CLINIC OPERATORY

- California Dental Chair
- 3 HPC Unit and 3 Way Syringe
- Telescoping Assistants Package
- Diamond Dental Light
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY



Add a Cuspidor



America Operator



AMERICA OPERATORY

- AM3 Dental Chair
- 3 HPC Unit and 3 Way Syringe
- Telescoping Assistants Package
- Diamond Dental Light
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY



10 Year
Warranty



Add a Cuspidor



Universal Operator

10 Year
Warranty

UNIVERSAL OPERATORY

- AM3 Dental Chair
- Radius 3 HPC Unit and 3 Way Syringe
- Rear Radius Assistants Package
- Radius Diamond Dental Light
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY



COMPLETE LIST OF ACCESSORIES ON PAGE 17



Independence Operator

10 Year
Warranty

INDEPENDENCE OPERATORY

- Infinity Dental Chair
- Cabinet Mounted 3 HPC Unit,
3 Way Syringe W/ Assistants Package
- Height Adjustable Work Surface
- Radius Diamond Dental Light
- Radius Monitor Bracket
- 800 Lift Capacity
- Over 150 Color Choices
- 10 YEAR WARRANTY

SEE THE NEW DANSEREAU
INDEPENDENCE CATALOG
FOR A COMPLETE GUIDE TO
THE FULL LINE OF REAR AND SIDE
DELIVERY SYSTEMS



COMPLETE LIST OF ACCESSORIES ON PAGE 17

OPERATOR SEATING

SIGNATURE DENTISTS STOOL



Ergonomic Controls
Aseptic Cushions
Aluminum Star
Chrome Cylinder

SIGNATURE ASSISTANTS STOOL



Adjustable Foot Ring
Aseptic Cushions
Aluminum Star
Chrome Cylinder

BASIC DENTISTS STOOL



Single Lever Lift
Aseptic Cushions
Aluminum Star
Chrome Cylinder

Ask about our Stool Package Discounts

CALIFORNIA DENTISTS STOOL



Single Lever Lift
Aseptic Cushions
Poly Resin Star
Black Cylinder
Non Marr Casters

CALIFORNIA ASSISTANTS STOOL



Adjustable Foot Ring
Aseptic Cushions
Poly Resin Star
Black Cylinder
Non Marr Casters

OPERATORS DENTISTS STOOL



Adjustable Foot Ring
Aseptic Cushions
Poly Resin Star
Black Cylinder
Non Marr Casters

12'O'CLOCK DELIVERY SYSTEMS

INDEPENDENCE CABINET / WALL SYSTEM



- 3 HPC Automatic Unit with 3 Way Syringe
- Complete Assistants package with HVE and Saliva Ejector
- Height Adjustable Work-surface
- Wall or Cabinet Installation
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY



Configures for Right or Left Use

UNIVERSAL CART SYSTEM



- 3 HPC Automatic Unit with 3 Way Syringe
- Complete Assistants package with HVE and Saliva Ejector
- Height Adjustable U Cart with Work-surface and EZ Roll Casters
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY



Configures for Right or Left Use

INFINITY CART SYSTEM



- 3 HPC Automatic Unit with 3 Way Syringe
- Complete Assistants package with HVE and Saliva Ejector
- Height Adjustable U Cart with EZ Roll Casters
- Unit Foot Control
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY



INFINITY DELIVERY SYSTEMS

INFINITY 3 DELIVERY SYSTEM

- 3 Handpiece Automatic System
- *Dentists* 3 Way Syringe Included
- Recessed Tray
- Option View Box (Add \$)
- Straight Soft Tubing
- Flex Arm with Air Lock
- Telescoping Assistants Arm
- *Assistants* 3 Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- Mounting Post and Complete Junction Box
- 10 YEAR WARRANTY



INFINITY 3 DELIVERY SYSTEM



- 3 Handpiece Automatic System
- *Dentists* 3 Way Syringe Included
- Recessed Tray
- Option View Box (Add \$)
- Straight Soft Tubing
- Flex Arm with Air Lock
- Cuspidor with Timed Rinse Bowl
- Manual Cupfiller
- Telescoping Assistants Arm
- *Assistants* 3 Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- Mounting Post and Complete Junction Box
- 10 YEAR WARRANTY

INFINITY 2 DELIVERY SYSTEM



- 3 Handpiece Automatic System
- 3 Way Syringe Included
- Recessed Tray
- Straight Soft Tubing
- Flex Arm with Air Lock
- Optional View Box (Add \$)
- 10 YEAR WARRANTY

INFINITY DELIVERY SYSTEM WILL INSTALL
ON A CHAIR, WALL OR CABINET.
ADVISE WHEN ORDERING

INFINITY ASSISTANTS SYSTEMS

INFINITY 1# ASSISTANTS SYSTEM

- Full Telescoping Assistants Arm
- Post Mounting Bracket
- Assistants 3 Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY



CUSPIDOR SYSTEM



- Cuspidor - Ceramic Bowl
- Short Telescoping Assistants Arm
- Timed Rinse Feature
- Manual Cupfiller
- Assistants 3 Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

INFINITY 2# TELESCOPING ARM SYSTEM



- Short Telescoping Assistants Arm
- Post Mounting Bracket
- Assistants 3 Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

10 Year
Warranty

INDEPENDENCE ASSISTANTS SYSTEM

- Radius Assistants Package
- Height Adjustable Laminated Table
- Assistants 3 Way Syringe
- 21 HVE Valve with Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY



INFINITY DELIVERY SYSTEM WILL INSTALL ON ANY CHAIR, WALL OR CABINET.
ADVISE WHEN ORDERING

INFINITY CONTROL SYSTEMS

INFINITY SLIDE OUT CONTROL UNIT W ASST. PACKAGE

INCLUDES UTILITY CENTER



- 3 Handpiece Automatic System
- *Dentists* 3 Way Syringe Included
- Straight Soft Tubing
- Complete Junction Box

10 Year
Warranty



INFINITY CONTROL UNIT W ASST. PACKAGE

INCLUDES UTILITY CENTER



- 3 Handpiece Automatic System
- *Dentists* 3 Way Syringe Included
- Straight Soft Tubing
- Complete Junction Box
- 10 YEAR WARRANTY

10 Year
Warranty



INFINITY FACE PLATE DELIVERY SYSTEM

INCLUDES UTILITY CENTER

- 2 Handpiece Automatic System
- *Dentists* 3 Way Syringe Included
- Straight Soft Tubing
- 1 HVE Valves with 1 Saliva Ejector
- Complete Junction Box
- 10 YEAR WARRANTY

10 Year
Warranty



INFINITY DELIVERY SYSTEM WILL INSTALL ON A CHAIR, WALL OR CABINET.
ADVISE WHEN ORDERING

EQUIPMENT ACCESSORIES

CHAIR FOOT CONTROL



CHAIR SWIVEL



ARTICULATING HEADREST NEW PUSH BUTTON



SELF CONTAINED WATER SYSTEM

FIBER OPTIC SYSTEM

4 Hole Connection
Traditional Fiber Optic Light
Source with the Bulb in the
Connector



6 Pin Connection
European Type Fiber Light
Source with the Bulb in the
Handpiece



UTILITY TRAY



COMPUTER MONITOR BRACKET



WET DRY FOOT CONTROL



WATER SYSTEM QUICK DISCONNECT



INFINITY PIEZO SCALE - BUILT IN



AIR COMPRESSOR SYSTEMS


CLEAN FLOW AIR COMPRESSOR

- 1 HP Oil Less Motor
- 120 volt or 240 volt
- Regenerating Drying System
- Automatic Purge
- 99.9% Coalescing Filter
- Installation Hardware Included
- Quiet Operation 62db
- Limited 5 Year Warranty



CLEAN FLOW AIR COMPRESSOR

- 2 HP Dual Oil Less Motor
- 240 volt
- Regenerating Drying System
- Automatic Purge
- 99.9% Coalescing Filter
- Installation Hardware Included
- Quiet Operation 62db
- Limited 5 Year Warranty



UTILITY ROOM ACCESSORIES

Low Voltage
Control
Panel



Water
Solenoid
Filter



Low Voltage
Control
Box



LARGE SIZE COMPRESSORS

CLEAN FLOW 3 HP (TRIPLE HEAD) COMPRESSOR WITH DRYER

CLEAN FLOW 4 HP (QUAD HEAD) COMPRESSOR WITH DRYER

VACUUM SYSTEMS

PRECISION VACUUM SYSTEM

- 1 HP Solid State Motor
- 120 volt or 240 volt
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add a Air / Water Separator
- Limited 5 Year Warranty



Note: 3 and 4 Horsepower Systems Available

PRECISION VACUUM SYSTEM

- 2 HP Dual Solid State Motors
- 240 volt only
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add an Air Water / Separator
- Limited 5 Year Warranty



SUPER FLOW 10

- 2 Horsepower
- Noise Level 40db
- Voltage - 240 volt
- Amps-6.4
- Vacuum in / HG - 7
- Stainless Steel Tank
- Tank Wash Included
- Versatile Installation
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty



ROADRUNNER - HYBRID

- 2 Horsepower
- Variable Speed Drive
- Adjusts to Vacuum Needs
- Noise Level 40db
- Voltage - 240 volt / Amps-6.4
- Vacuum in / HG - 7
- Stainless Steel Tank
- Wash Down System Included
- Reduced Electrical Costs
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty



Real Time Vacuum
Adjustments based upon
Vacuum Demands

DENTAL LIGHTS

INFINITY LED DENTAL LIGHT

INSTALLS ON ANY DENTAL CHAIR

- Light Head with 3 Axis Movement
- LED Light Source
- Adjustable Intensity
- Color Temperature: 4800 kelvin
- Color Rendering Index 92
- Red Coloring Index 60
- Focal Shape Oval
- Counter Balance, Drift Free Arm



3 Axis Rotation



Dimension Specifications

Length of Horizontal Arm: 34"

Length of Flex Arm: 26"

Flex Arm Vertical Movement: 31"

Includes Transformer

DIAMOND DENTAL LIGHT 2014

INSTALLS ON ANY DENTAL CHAIR

- Light Head with 2 Axis Movement
- Xenon Quartz Halogen Bulb
- Dual Intensity - High / Low Illumination
- Color Temperature: 4200 to 5000 kelvin
- Counter Balance, Drift Free Arm
- Ergonomic Handles
- Simple Lever Activation
- Glass Reflector with Shadow Free Illumination



Dimension Specifications

Length of Horizontal Arm: 24"

Length of Flex Arm: 26"

Flex Arm Vertical Movement: 31"

Includes Transformer

STERILIZERS / ULTRASONICS

TUTTNAUER 1730 STERILIZER

- Chamber Size - 7" x 13"
- Cycle Time 11 Minutes (Hot Start)
- Auto Shut Off on Both Cycles
- 120 volt - 11amp
- Easy Drain for Water Reservoir
- Tray Capacity - 2 Half
- 1 Year Warranty Parts and Labor



TUTTNAUER 2340M STERILIZER

- Chamber Size - 9" X 18"
- Cycle Time 13 Minutes (Hot Start)
- Auto Shut Off on Both Cycles
- 120 volt - 12amp
- Easy Drain for Water Reservoir
- Tray Capacity - 2 Full & 2 Half
- 1 Year Warranty Parts and Labor



TUTTNAUER EZ 9 STERILIZER

- Chamber Size - 9" X 18"
- Cycle Time 13 Minutes (Hot Start)
- Auto Drying Cycle
- Low Water Sensor
- Pre-programmed Sterilization Cycles
- 120 volt - 12amp
- Easy Drain for Water Reservoir
- Tray Capacity - 2 Full & 2 Half



TUTTNAUER EZ 10 STERILIZER

- Chamber Size - 10" X 19"
- Cycle Time 14 Minutes (Hot Start)
- Auto Drying Cycle
- Low Water Sensor
- Pre-programmed Sterilization Cycles
- 120 volt - 12amp
- Easy Drain for Water Reservoir
- Tray Capacity - 2 Full & 2 Half



TUTTNAUER DISTILLER

- 1 Gallon



TUTTNAUER CSU1 ULTRASONIC

- 1 Gallon
- 2 Cassettes



TUTTNAUER CSU3 ULTRASONIC

- 3 Gallon
- 3 Cassettes



X-RAY MACHINES

BELMONT PHOT-XII INTRAORAL X-RAY MACHINE

CALL FOR SPECIAL PRICING

- Wall or Cabinet Mounted
- 70kvp - 8mA
- Warranty 2 years
- Extension 75.5"
- Automatic Exposure Times
- Digital or Standard X - rays




BELMONT 097 BELRAY INTRAORAL X-RAY MACHINE

CALL FOR SPECIAL PRICING

- Wall or Cabinet Mounted
- 70kvp - 8mA
- Warranty 2 years
- Extension 69"
- 12" Long Cone
- Digital or Standard X - rays



DENTAL OFFICE EQUIPMENT PACKAGE



Typical Equipment List:

Low Voltage Control Panel / Low Voltage Transformer / DHP 1 HP Oil - Less Air Compressor
 DHP 1 HP Vacuum Pump / DHP Air / Water Separator / 2# - America Operatories with Stools
 Intra Oral X ray / Velopex X - ray Processor / Tuttnauer 1730 Sterilizer.

Delivery and Installation optional.

X - ray Processor Not Included*

**NOTE: ALL THIS EQUIPMENT IS NEW
AND FULLY WARRANTIED**

RXDigital

RX Soluciones de Imagenología

RX SOLUCIONES DE IMAGENOLOGÍA S.A. de C.V.

Av. Moctezuma 3515-1B, Col. Ciudad del Sol. Zapopan, Jalisco. MÉXICO

 VISITE Nuestro Sitio Web
 y vea nuestra Línea de
 Productos Completa


TEL (33) 3880-1262 CEL (33) 3499-7095 / 1300-5813

contacto@rxdigitalsoluciones.com

ventas@rxdigitalsoluciones.com

www.rxdigitalsoluciones.com



RxDigitalDental

- Design your own Dental Office Package: You can add or eliminate any products from the package.
- Office Design is also available. Simple Designs to Complete Blueprints.

DANSEREAU DENTAL

AMERICA'S DENTAL EQUIPMENT COMPANY

Dansereau has the most comprehensive product warranty programs available to the dentist today. We offer a 10 Year Warranty on our Dental Chairs and Delivery Systems with complete features listed below:

5 Years 100% / 6-10 Years Prorated Annually

Dental Chair Electric Motors, Chair Base Frame and Chair Frame.

3 Years 100% / 4-10 Years Prorated Annually

Delivery System Internal Tubing (against leaks), Master Valve Regulators, Handpiece Control Blocks (Not Gaskets), Light Posts, Automatic Handpiece Holders, Delivery System Water and Air Flow Control Valves, Stool Cylinders, Manual Chair Switches, Delivery System Flex Arm, Delivery System Telescoping Arm, Cuspidor Frame.

1 Year 100%

Diamond Dental Light, Chair Foot Control, Stool Ring, Stool Seat, Stool Casters, Circuit Board and Prepositioning Switches.

6 Months 100%

Base and Unit Cover parts, 3-way Syringes, Delivery System Foot Control and Fiber Optics Systems.

30 Days 100%

Dental Light Bulbs, Fiber Optic Light Bulbs, X-ray View Box Light Bulbs.

Broken or User Damaged Items are Not Covered under Warranty.

Minimum Maintenance is required in order to guarantee full use of warranties

All dental delivery systems require clean dry air to allow the proper lifespan to the micro valves within the dental delivery system. If your air compressor does not have an adequate drying and filtration system the warranties of the micro valves, on/off toggle valves, relays, regulators, master valves, foot controls, 3 ways syringes and automatic blocks may be void.

GUARANTEE

Dansereau will guarantee against defects in material and workmanship. Defective parts will be repaired or replaced: no claim for labor or consequential damages will be allowed. All replacement parts claimed defective will be filled in the usual manner, such invoice being subject to adjustment after the alleged defective material has been returned and inspected by the manufacturer. The warranty does not cover damage resulting from improper installation, maintenance, accident or misuse.

Terms: Prepaid / Financing is Available

Freight: All Shipments are FOB Corona, CA. unless otherwise specified by purchaser on the order. Dansereau is not responsible for any freight unless otherwise specified in writing.

Returned Product: All orders are Non Cancelable, No parts can be returned without a factory authorization. All parts are subject to 25% re-stocking fee and must be returned Prepaid Shipping.